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
UCMSC, α -MEM

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
The total cell number, compared with 20% Defined FBS

5% UltraGRO™ 3X↑

Doubling Time 25 hours

Gwowe, 2013


BMSC, α -MEM

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
The total cell number, compared with 20% Defined FBS

5% UltraGRO™ 3X↑

Doubling Time 28 hours

GwoWei, 2013


ADMSC, α -MEM

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
The total cell number, compared with 20% FBS

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




▪ *Fibroblast growth factor-2 (FGF2), platelet-derived growth factor (PDGF), epidermal growth factor (EGF), and transforming growth factor (TGF)- β are the growth factors which are important for MSCs proliferation.*

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


UltraGRO®-PURE GI, GMP grade


Fig. 2: Growth factors retained comparable cytokine levels after receiving gamma irradiation.

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
Fig.3: Human adipose tissue derived MSCs retain tri-lineage differentiation capability after cultured in UltraGRO™-PURE GI supplemented medium for three passages.


Fig.4: MSCs retained immunomodulation potency

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Cell Reports Medicine


Article

iPSC-derived hypoimmunogenic tissue resident memory T cells mediate robust anti-tumor activity against cervical cancer

Yoshiki Furukawa,^{1,8} Midori Ishii,^{1,8} Jun Ando,^{1,2} Kazuya Ikeda,³ Kyomi J. Igarashi,⁴ Shintaro Kinoshita,¹ Yoko Azusawa,² Tokuko Toyota,¹ Tadahiro Honda,¹ Mahito Nakanishi,⁵ Koichi Ohshima,⁶ Ayako Masuda,⁷ Emiko Yoshida,⁷ Mari Kitade,⁷ Matthew Porteus,³ Yasuhisa Terao,⁷ Hiromitsu Nakauchi,^{4,*} and Miki Ando^{1,9,*}

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Article

Interferon- α -Induced Dendritic Cells Generated with Human Platelet Lysate Exhibit Elevated Antigen Presenting Ability to Cytotoxic T Lymphocytes

Ippei Date ¹, Terutsugu Koya ^{1,2}, Takuya Sakamoto ^{1,2}, Misa Togi ^{1,2}, Haruhiko Kawaguchi ¹, Asuka Watanabe ¹, Tomohisa Kato Jr. ³ and Shigetaka Shimodaira ^{1,2,*}

Vaccines 2021,9(1), 10

<https://dx.doi.org/10.3390/vaccines9010010>

Article

Clinical Trial on the Safety and Tolerability of Personalized Cancer Vaccines Using Human Platelet Lysate-Induced Antigen-Presenting Cells

Terutsugu Koya ^{1,2}, Kenichi Yoshida ², Misa Togi ^{1,2}, Yo Niida ³, Sumihito Togi ³, Hiroki Ura ³, Shuichi Mizuta ⁴, Tomohisa Kato, Jr. ⁵, Sohsuke Yamada ⁶, Takeo Shibata ⁷, Yi-Chang Liu ^{8,9}, Shyng-Shiou Yuan ¹⁰, Deng-Chyang Wu ¹¹, Hirohito Kobayashi ¹², Taiju Utsugisawa ¹³, Hitoshi Kanno ¹³ and Shigetaka Shimodaira ^{1,2,5,13,*}

Cancers 2023, 15(14), 3627

<https://doi.org/10.3390/cancers15143627>